



Leica EM VCT100

Vacuum Cryo Transfer

From Preparation to Analysis

Living up to Life

Leica

MICROSYSTEMS

Versatile Specimen Transfer

Transferring specimens to the chamber of an analysis system is a very critical step for most preparation methods. It is essential to protect the samples from contamination. The Leica EM VCT100 concept was established to cross-link preparation units with various analysis systems via a transfer shuttle connected to a load-lock. Samples are transferred in a well-defined environment, e.g. protective gas or high vacuum conditions from a preparation instrument to the analysis unit. In addition, the specimen can be kept at low temperature for cryo techniques. The load-lock can mount to any vacuum operated SEM, FIB, SIMS, AFM or XPS chambers.

The sample holders used during specimen preparation and transfer are compatible with the custom-made stage of the analysis instrument. Due to the modular concept, the stage can be installed in any (ultra) high vacuum system or glove box and used for both room temperature and cryo applications. The flexible system design enables adjustments to any specific application.



EM VCT100 shuttle with specimen and holder



EM VCT100 cryo preparation workstation

Setup Advantages and Key Features

- Contamination-free transfer between preparation and analysis unit
- Unique shuttle system for the vacuum cryo transfer of specimens
- Special shuttle and load-lock design maintains resolution of the SEM, vibration free
- Suitable for either room or cryo transfer
- Versatile protective gas transfer such as argon to prevent oxidation
- Preparation and analysis can be performed at different locations
- Space saving design offers minimal interference with the analysis system
- Sample preparation and analysis can be performed independently without interruption of either process
- Possible to repeat preparation with the same specimen
- Adaptable to more than one SEM
- Preparation units can be linked to several Leica EM VCT100 adapted analysis units such as the EM SCD500, EM MED020 or the EM BAF060

Leica EM VCT100 control unit enables integrated safety communication for specimen transfer through the load-lock system of the preparation and analysis unit.



Freeze etch/fracture unit with the Leica EM VCT100 shuttle attached. From preparation to imaging, the EM VCT100 is a complete cryo transfer system. Its components include among others a cryo preparation workstation, transfer shuttle, docking station, a controlled cold stage (-150°C to +60°C), LN₂ Dewar and operating panel with touch screen. For docking to a SEM or different analysis instrument such as FIB, SIMS, XPS, etc., a cryo adaption kit with cryo stage is included.

The Leica EM VCT100 for Cryo SEM

Cryo Electron Microscopy (Cryo SEM) is the process of observing non-stained and untreated vitreous samples. Cryo sample preparation is a quick sample preparation technique that has none of the artefacts associated with chemical fixation. It requires no toxic reagents and reduces beam damage due to low temperature

during imaging. With the right experimental setup the entire processing of the specimen from excision to immobilization (freezing), surface preparation, coating, and imaging can be done in less than one hour.

STEP 1: Sample Freezing

A typical preparation procedure starts with cryo fixation of the specimen using techniques like high pressure freezing with the Leica EM HPM100 or the Leica EM PACT2, jet freezing with the Leica EM JFD 030 or plunge freezing with the Leica EM CPC or Leica EM GP, depending on the nature of the specimen and application desired.



Leica EM HPM100



Leica EM PACT2



Leica EM GP




Leica EM CPC

STEP 2: Sample Loading

The specimen is loaded under LN₂ onto a specimen holder using the loading chamber of the Leica EM VCT100 loading box and then retracted into the pre-cooled shuttle of the EM VCT100. Since the shuttle and loading chamber are pre-cooled with LN₂, the specimen is transferred to the preparation unit contamination-free without any exposure to air.

There are a variety of specimen holders to choose from, e.g. for freeze fracture, double replica, freeze drying, etc.



STEP 3: Sample Transfer for Fracturing, Etching and Coating

The specimen and holder are subsequently transferred under vacuum and low temperature to one of the many Leica preparation systems using the Leica EM VCT100 shuttle. The specimen can then be freeze fractured, freeze etched, freeze dried and/or coated for follow-on analysis in the cryo SEM.

Three preparation instruments are available for cryo SEM sample preparation: the Leica EM SCD500, the EM MED020 and the EM BAF060.

Leica EM SCD500

The Leica EM SCD500 offers many conversion options in one single unit, all of which can be easily adapted to a variety of applications:

- high vacuum sputtering
- single and multiple carbon thread evaporation
- thermal resistance evaporation
- carbon rod evaporation
- etching and glow discharge
- cryo coating, double replica, freeze drying, freeze etching and vacuum cryo transfer with the Leica EM VCT100



Leica EM MED020

The modular high vacuum system can be fitted with different, optional attachments (i.e., modules) for the following sample preparations:

- single and triple high vacuum sputter coating
- single and multiple carbon thread evaporation
- thermal resistance evaporation
- carbon rod evaporation
- e-beam evaporation
- etching and glow discharge
- cryo preparation for freeze drying, freeze fracturing, double replica, freeze etching, cryo coating and vacuum cryo transfer with the EM VCT100



Leica EM BAF060


The Leica EM BAF060 is a fully automatic high-end preparation unit for:

- freeze fracturing
- freeze etching
- freeze drying
- double replica (mirror fracturing)
- high resolution carbon/metal mix coatings for TEM/SEM analysis
- specimen replication by electron beam evaporation
- double layer coated specimens for cryo SEM analysis
- cryo coating and vacuum cryo transfer with the Leica EM VCT100




STEP 4: Sample Transfer for Analysis

After sample preparation, the specimen is transferred still frozen and under vacuum via the Leica EM VCT100 shuttle onto the cryo stage within the SEM.



EM VCT100 docking station and shuttle with specimen and holder connected to the cryo stage.

The custom-made cryo stage slides onto the standard SEM stage and is connected via copper bands to the LN₂ Dewar. The specimen holder with the specimen is transferred onto the cryo stage using the EM VCT100 shuttle. The cryo stage can be used for any specimen with the advantage of improving the vacuum in the chamber and reducing beam damage. SEM resolution remains unchanged. The stage can also be used for room temperature analysis with the cryo stage in position.



Custom-made cryo stage

During SEM analysis the Leica EM VCT100 shuttle is detached from the SEM so it has no effect on the stability of the microscope.




The EM VCT100 is compatible with all current SEMs.

Specimen Holders

A variety of specimen holders can be supplied with the system. The selection depends on the cryo fixation technique and the specimen.

① **Freeze Fracture holder**

with retaining spring for three $\varnothing 3 \times 0.8$ mm gold specimen carriers.



② **Freeze Fracture holder**

with retaining spring for three $\varnothing 3 \times 4.5 \times 0.8$ mm copper specimen carriers.



④ **Double Replica holder**

with spring load for two $\varnothing 4.6 \times 0.6$ mm gold specimen carrier sandwiches.



⑤ **Double Replica holder**

with spring load for four $\varnothing 3 \times 4.5 \times 0.6$ mm copper specimen carrier sandwiches.



⑥ **Freeze Drying holder**

with two magnetic strips for three $\varnothing 3$ mm grids.

⑦ **Freeze Drying holder**

with retaining spring for four $\varnothing 3$ mm grids.

⑧ **Holder for SEM stubs**

with fastening thread (M3) for screw-on SEM stubs.


⑨ **Blank holder**

The specimen is adhered directly onto the holder.




Applications

Cryo SEM is an imaging technique for a wide range of applications in life and material sciences. Cryo SEM is not exclusively used to access the ultra-structure of an object. This technique can also be used for the analysis and localization of different components (e.g. elements).




High pressure frozen, freeze fractured cement suspension. The sample was fractured at -115°C in a BAF060, etched at -105°C for 4 minutes and coated by electron beam evaporation with 3 nm of Pt/C. P...Cement particle, E...Ettringites. Electron Microscopy ETH Zurich (EMEZ).

Courtesy of L. Holzer, EMPA, Swiss Federal Laboratories for Materials Testing and Research, 3D-Mat group, Duebendorf, Switzerland.




High pressure frozen, freeze fractured liposome suspension. The suspension was fractured at -150°C in a BAF060 and double layer coated by electron beam evaporation with 3 nm of Pt/C at 45° and subsequently with 6 nm of carbon. Specimen was imaged using backscattered electrons. Electron Microscopy ETH Zurich (EMEZ).




High pressure frozen, freeze fractured mouse intestine biopsy. The sample was fractured at -115°C in a BAF060, etched at -105°C for 5 minutes and coated by electron beam evaporation with 3 nm of Pt/C. M...Mitochondria, N...Nucleus, NP...Nucleopores, G...Golgi, PM...Plasma membrane. Electron Microscopy ETH Zurich (EMEZ).

Courtesy of René Fischer, Laboratory of Organic Chemistry, ETH Zurich, Switzerland.




High pressure frozen, freeze fractured suspension of Vero cells. The sample was fractured at -115°C in a BAF060 and immediately coated by electron beam evaporation with 3 nm of Pt/C. M...Mitochondria, N...Nucleus, NP...Nucleopores, G...Golgi, PM...Plasma membrane. Electron Microscopy ETH Zurich (EMEZ).

Courtesy of Peter Wild, Institute of Veterinary Anatomy, University of Zurich, Switzerland.



High pressure frozen (HPM010), EM VCT100 transfer to the EM BAF060 for freeze-fracture/freezing-etching and cryo-coating (Pt/C, 3 nm) using the electron beam gun and rotating specimen holder. EM VCT100 transfer to the cryo-FEGSEM (JEOL 7401F). *Euglena gracilis* Klebs CCAP 1224/5Z.

Courtesy of Dr. Roland Fleck, NIBSC, Potters Bar, UK.



High pressure frozen (HPM010), EM VCT100 transfer to the EM BAF060 for freeze-fracture/freezing-etching and cryo-coating (Pt/C, 3 nm) using the electron beam gun and rotating specimen holder. EM VCT100 transfer to the cryo-FEGSEM (JEOL 7401F). Pennate diatom from a mixed culture of the protist Euplotes.

Courtesy of Dr. Roland Fleck, NIBSC, Potters Bar, UK.

“With the user, for the user”

Leica Microsystems

Leica Microsystems operates globally in four divisions, where we rank with the market leaders.

• Life Science Division

The Leica Microsystems Life Science Division supports the imaging needs of the scientific community with advanced innovation and technical expertise for the visualization, measurement, and analysis of microstructures. Our strong focus on understanding scientific applications puts Leica Microsystems' customers at the leading edge of science.

• Industry Division

The Leica Microsystems Industry Division's focus is to support customers' pursuit of the highest quality end result. Leica Microsystems provide the best and most innovative imaging systems to see, measure, and analyze the microstructures in routine and research industrial applications, materials science, quality control, forensic science investigation, and educational applications.

• Biosystems Division

The Leica Microsystems Biosystems Division brings histopathology labs and researchers the highest-quality, most comprehensive product range. From patient to pathologist, the range includes the ideal product for each histology step and high-productivity workflow solutions for the entire lab. With complete histology systems featuring innovative automation and Novocastra™ reagents, Leica Microsystems creates better patient care through rapid turnaround, diagnostic confidence, and close customer collaboration.

• Surgical Division

The Leica Microsystems Surgical Division's focus is to partner with and support surgeons and their care of patients with the highest-quality, most innovative surgical microscope technology today and into the future.

The statement by Ernst Leitz in 1907, “with the user, for the user,” describes the fruitful collaboration with end users and driving force of innovation at Leica Microsystems. We have developed five brand values to live up to this tradition: Pioneering, High-end Quality, Team Spirit, Dedication to Science, and Continuous Improvement. For us, living up to these values means: **Living up to Life.**

Active worldwide

Australia:	North Ryde	Tel. +61 2 8870 3500	Fax +61 2 9878 1055
Austria:	Vienna	Tel. +43 1 486 80 50 0	Fax +43 1 486 80 50 30
Belgium:	Groot Bijgaarden	Tel. +32 2 790 98 50	Fax +32 2 790 98 68
Canada:	Richmond Hill/Ontario	Tel. +1 905 762 2000	Fax +1 905 762 8937
Denmark:	Herlev	Tel. +45 4454 0101	Fax +45 4454 0111
France:	Nanterre Cedex	Tel. +33 811 000 664	Fax +33 1 56 05 23 23
Germany:	Wetzlar	Tel. +49 64 41 29 40 00	Fax +49 64 41 29 41 55
Italy:	Milan	Tel. +39 02 574 861	Fax +39 02 574 03392
Japan:	Tokyo	Tel. +81 3 5421 2800	Fax +81 3 5421 2896
Korea:	Seoul	Tel. +82 2 514 65 43	Fax +82 2 514 65 48
Netherlands:	Rijswijk	Tel. +31 70 4132 100	Fax +31 70 4132 109
People's Rep. of China:	Hong Kong	Tel. +852 2564 6699	Fax +852 2564 4163
Portugal:	Lisbon	Tel. +351 21 388 9112	Fax +351 21 385 4668
Singapore		Tel. +65 6779 7823	Fax +65 6773 0628
Spain:	Barcelona	Tel. +34 93 494 95 30	Fax +34 93 494 95 32
Sweden:	Kista	Tel. +46 8 625 45 45	Fax +46 8 625 45 10
Switzerland:	Heerbrugg	Tel. +41 71 726 34 34	Fax +41 71 726 34 44
United Kingdom:	Milton Keynes	Tel. +44 1908 246 246	Fax +44 1908 609 992
USA:	Bannockburn/Illinois	Tel. +1 847 405 0123	Fax +1 847 405 0164

and representatives in more than 100 countries